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
9. DOCUMENTOS DE REFERENCIA

- Bates, D. M. and Watts, D. G.: Nonlinear Regression : Iterative Estimation and Linear Approximations, in Nonlinear Regression Analysis and Its Applications, JOHN WILEY & SONS, Ontario, Canada., 1988.
- Dickson, A. G.: An exact definition of total alkalinity and a procedure for the estimation of alkalinity and total inorganic carbon from titration data, *Deep Sea Res. Part A, Oceanogr. Res. Pap.*, 28(6), 609–623, doi:10.1016/0198-0149(81)90121-7, 1981.
- Dickson, A. G.: The development of the alkalinity concept in marine chemistry, *Mar. Chem.*, 40(1–2), 49–63, doi:10.1016/0304-4203(92)90047-E, 1992.
- Dickson, A. G., Sabine, C. L. and Christian, J. R.: Guide to best practices for ocean CO₂ measurements, *PICES Spec. Publ.* 3, 3(8), 191, doi:10.1159/000331784, 2007.
- Ellison, S. L. . and Williams, A.: *Eurachem/CITAC guide: Quantifying Uncertainty in Analytical Measurement*, Third edit., 2012.
- Gazeau, F., Alliouane, S. and Comeau, S.: Total alkalinity Protocol, *Observatoire Oceanologique Villefranche sur mer. Francia.*, 2014.
- Gran, G.: Determination of the equivalence point in potentiometric titrations. Part II, *Analyst*, 77(920), 661–670, doi:10.1039/AN9527700661, 1952.
- Hernández-Ayon, J. M., Zirino, A., Dickson, A. G., Camiro-Vargas, T. and Valenzuela-Espinoza, E.: Estimating the contribution of organic bases from microalgae to the titration alkalinity in coastal seawaters, *Limnol. Oceanogr. Methods*, 5(JUL), 225–232, doi:10.4319/lom.2007.5.225, 2007.
- Huang, W. J., Wang, Y. and Cai, W. J.: Assessment of sample storage techniques for total alkalinity and dissolved inorganic carbon in seawater, *Limnol. Oceanogr. Methods*, 10, 711–717, doi:10.4319/lom.2012.10.711, 2012.
- Lavigne, H., Orr, J., Gentili, B., Hagens, M., Hofmann, A. and Rae, J.: Package ‘ seacarb ’: Seawater Carbonate Chemistry. [online] Available from: <https://cran.r-project.org/package=seacarb>, 2019.
- Metrohm: 888 Titrand. User manual, 8.888.8004. [online] Available from: www.methom.com, 2016a.
- Metrohm: TiamoTM. Tutorial, 8.101.8003. [online] Available from: www.methom.com, 2016b.
- Pimenta, A. . and Grear, J. .: Guidelines for measuring changes in seawater pH and associated carbonate chemistry in coastal environments of the Eastern United States, p. 56, Office of Research and Development, National Health and Environmental Effects Research Laboratory., 2018.
- Wolf-Gladrow, D. A., Zeebe, R. E., Klaas, C., Körtzinger, A. and Dickson, A. G.: Total alkalinity: The explicit conservative expression and its application to biogeochemical processes, *Mar. Chem.*, 106(1-2 SPEC. ISS.), 287–300, doi:10.1016/j.marchem.2007.01.006, 2007.

10. ANEXOS

10.1. ANEXO 1: Configuración del método “Verification with NBS Buffers”

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The screenshot displays the 'tiarno 2.5 - Workplace' software interface. On the left, a vertical toolbar contains icons for Workplace, Database, Method, Configuration, and Manual. A central workflow diagram shows the following steps: START (Main track), CAL LOOP pH (Calibration loop pH), CAL MEAS pH (CAL MEAS pH), REPORT (Report), DATABASE (Transfer to database), and END (Main track). Three configuration dialog boxes for the 'START - Main track' command are overlaid on the right side of the screen.


The first dialog box, titled 'START - Main track', shows the 'General' tab with the following settings:

- Command name: Main track
- Workplace view: Current view, View
- Track view for live window: Live display 1: Main track, Live display 2: Main track
- Statistics: Statistics
- Number of single determinations: not defined
- Conditioning: Automatic conditioning


The second dialog box, also titled 'START - Main track', shows the 'Application note' tab, which is currently empty.

The third dialog box, titled 'START - Main track', shows the 'Method variables' tab, which contains a table with the following columns: Name, Type, Assignment, Fixed value, Comment, and Monitoring. The table is currently empty.

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
The screenshot displays the 'tiamo 2.5 - Workplace' software interface. On the left, a vertical toolbar contains icons for Workplace, Database, Method, Configuration, and Manual. The main workspace shows a workflow diagram with the following steps: START (Main track), CAL LOOP pH (Calibration loop pH), CAL MEAS pH (CAL MEAS pH), REPORT (Report), DATABASE (Transfer to database), and END (Main track). Two configuration windows for the 'CAL MEAS pH' command are open. The left window shows the 'General/Hardware' tab with the following settings: Command name: CAL MEAS pH; Device name: 888_1; Device type: 888 Titrande; Measuring input: 1; Sensor: pH electrode; Temperature measurement: automatic; Stirrer: 1; Stirring rate: 2; and a checked option for 'Switch off automatically'. The right window shows the 'Measuring parameters' tab with the following settings: Signal drift: 0.2 mV/min; Min. waiting time: 10 s; Max. waiting time: 332 s; Measuring interval: 2.0 s; and Temperature: 22 °C. Both windows have OK, Cancel, and Apply buttons at the bottom.

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
The screenshot displays the 'tiarno 2.5 - Workplace' software interface. On the left, a vertical toolbar contains icons for Workplace, Database, Method, Configuration, and Manual. The main workspace shows a workflow diagram with the following steps: START (Main track), CAL LOOP pH (Calibration loop pH), CAL MEAS pH (CAL MEAS pH), REPORT (Report), DATABASE (Transfer to database), and END (Main track). A 'REPORT - Report' dialog box is open, showing the following configuration options:

- Command name: Report
- Report template: calibration report
- Report output:
 - Printer: Default printer
 - PDF file: _exports\Calibration_electrode.pdf
 - Send e-mail: E-mail...

An arrow points to the PDF file path, with the text: "Indique ruta de acceso a directorio Calibration_Export".

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The screenshot displays the 'tiamo 2.5 - Workplace' interface. On the left, a vertical toolbar contains icons for 'Workplace', 'Database', 'Method', 'Configuration', and 'Manual'. The main workspace shows a workflow diagram for 'Method 'Verification electrode with NBS buffers' - Version 21'. The workflow consists of the following steps:

- START** (Main track)
- CAL LOOP pH** (Calibration loop pH)
- CAL MEAS pH** (CAL MEAS pH)
- REPORT** (Report)
- DATABASE** (Transfer to database)
- END** (Main track)

A dialog box titled 'DATABASE - Transfer to database' is open, showing a table with the following data:

| Database | |
|----------|-------------|
| 1 | Calibration |

Below the table are buttons for 'New', 'Properties', and 'Delete'. At the bottom of the dialog are 'OK', 'Cancel', and 'Apply' buttons. The status bar at the bottom left indicates 'Press F1 for help'.

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START - Main track

Command name: Main track

General Application note Method variables

Workplace view

Current view

View

Track view for live window

Live display 1: Main track

Live display 2: Main track

Statistics

Number of single determinations: not defined

Conditioning

Automatic conditioning

OK Cancel

START - Main track

Command name: Main track

General Application note Method variables

START - Main track

Command name: Main track

General Application note Method variables

| Name | Type | Assignment | Fixed value | Comment | Monitoring |
|------------|--------|------------|-------------|---------|--------------------------|
| 1 Salinity | Number | ID1 | | | <input type="checkbox"/> |

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The screenshot displays the 'tiamo 2.5 - Calibration electrode with TRIS buffer' software interface. On the left, a vertical toolbar includes icons for Workplace, Database, Method, Configuration, and Manual. The main workspace shows a workflow diagram with the following steps: START (Main track), MEAS U (MEAS U), CALC (Calcul mV, temp and pH), DATABASE (Transfer to database), EXPORT (Export Calibration_TRIS), and END (Main track).

Four configuration windows for the 'MEAS U' command are overlaid on the main interface:

- MEAS U - MEAS U (Top Left):** Shows the 'General/Hardware' tab. Fields include: Command name (MEAS U), Device name (888_1), Device type (888 Titrande), Measuring input (1), Sensor (pH electrode), Temperature measurement (automatic), and Stirrer (1, 2). A checkbox for 'Switch off automatically' is checked.
- MEAS U - MEAS U (Top Right):** Shows the 'Measuring parameters' tab. Fields include: Command name (MEAS U), Signal drift (0.2 mV/min), Min. waiting time (30 s), Max. waiting time (332 s), Measuring interval (1 s), Stop measured value (off mV), and Temperature (22 °C). The 'Measurement with drift control' radio button is selected.
- MEAS U - MEAS U (Bottom Left):** Shows the 'Evaluations' tab. Under 'Additional measured values', there are checkboxes for 'Additional calculated measured values' and 'Additional external measured values', both of which are currently unchecked.
- MEAS U - MEAS U (Bottom Right):** Shows the 'Evaluations' tab. Under 'Additional measured values', there are checkboxes for 'Fixed endpoint evaluation', 'Minimum evaluation', 'Maximum evaluation', and 'Break point evaluation', all of which are currently unchecked.

The Windows taskbar at the bottom shows the system clock as 00:22 on 26-May-20.

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tiamo 2.5 - Calibration electrode with TRIS buffer

File Edit View Insert Help

100 %

1 Workplace

1 Database

2 Method

Configuration

Manual

START
Main track

MEAS U
MEAS U

CALC
Calcul mV, temp and pH

DATABASE
Transfer to database

EXPORT
Export Calibration_TRIS

END
Main track

CALC - Calcul mV, temp and pH

Command name: Calcul mV, temp and pH

| Result name | Formula | Unit | Decimal places | Assignment | Statistics | Result monitoring |
|-------------|---------------------------------------|------|----------------|------------|-------------------------------------|--------------------------|
| 1 mv_TRIS | = 'MEAS U.EME' | | 2 | RS01 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| 2 Temp_TRIS | = 'MEAS U.ETE' | °C | 2 | RS02 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| 3 pH_TRIS | =(11911,08-18,2499* 'MV.Salinity' ... | | 2 | RS03 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |

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tiemo 2.5 - Calibration electrode with TRIS buffer

File Edit View Insert Help

100 %

1 Workplace

1 Database

2 Method

Configuration

Manual

START
Main track

MEAS U
MEAS U

CALC
Calcul mV, temp and pH

DATABASE
Transfer to database

EXPORT
Export Calibration_TRIS

END
Main track

CALC - Calcul mV, temp and pH

Command name: Calcul mV, temp and pH

| | Result name | Formula | Unit | Decimal places | Assignment | Statistics | Result monitoring |
|---|-------------|---------------------------------------|------|----------------|------------|-------------------------------------|--------------------------|
| 1 | mv_TRIS | = 'MEAS U.EME' | | 2 | RS01 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| 2 | Temp_TRIS | = 'MEAS U.ETE' | °C | 2 | RS02 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |
| 3 | pH_TRIS | =(11911,08-18,2499* 'MV.Salinity' ... | | 2 | RS03 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |

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The screenshot displays the 'tiamo 2.5' software interface for 'Calibration electrode with TRIS buffer'. The main window shows a workflow diagram with the following steps:

- START** (Main track)
- MEAS U** (MEAS U)
- CALC** (Calcul mV, temp and pH)
- DATABASE** (Transfer to database)
- EXPORT** (Export Calibration_TRIS)
- END** (Main track)

A dialog box titled 'DATABASE - Transfer to database' is open, showing the following configuration:

- Command name:
- Database list:

| | Database |
|---|-------------|
| 1 | Calibration |
- Buttons: New, Properties, Delete, OK, Cancel


The interface includes a sidebar with icons for Workplace, Database, Method, Configuration, and Manual. The status bar at the bottom indicates 'Press F1 for help'.

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The screenshot displays the 'tiamo 2.5 - Calibration electrode with TRIS buffer' software interface. On the left, a vertical toolbar contains icons for Workplace, Database, Method, Configuration, and Manual. The main workspace shows a vertical flowchart with the following steps: START (Main track), MEAS U (MEAS U), CALC (Calcul mV, temp and pH), DATABASE (Transfer to database), EXPORT (Export Calibration_TRIS), and END (Main track). An 'EXPORT - Export Calibration_TRIS' dialog box is open, showing 'Command name' as 'Export Calibration_TRIS' and 'Export template' as 'Calibration_export'. The dialog has 'OK' and 'Cancel' buttons at the bottom. The status bar at the bottom left reads 'Press F1 for help'.

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The screenshot displays the 'tiamo 2.5 - Alkalinity for 50 ml samples' application window. On the left, a vertical toolbar contains icons for Workplace, Database, Method, Configuration, and Manual. The main workspace shows a workflow diagram with the following steps: START (Main track), MEAS U (Initial pH-Temperature), CALC (Calcul Voltage pH 4), CALC (Calcul Voltage pH 3), MET U (Titration), and EXPORT (Export to csv). Three configuration dialog boxes for the 'STAR - Main track' command are open:

- General tab:** Command name is 'Main track'. Under 'Workplace view', 'Current view' is selected. Under 'Track view for live window', 'Live display 1' and 'Live display 2' are both set to 'Main track'. Under 'Statistics', 'Number of single determinations' is 'not defined'. Under 'Conditioning', 'Automatic conditioning' is checked.
- Application note tab:** This tab is currently empty.
- Method variables tab:** Contains a table with the following data:

| | Name | Type | Assignment | Fixed value | Comment | Monitoring |
|---|-----------------|--------|-------------|-------------|-------------|--------------------------|
| 1 | Sample name | Text | ID2 | | Sample name | <input type="checkbox"/> |
| 2 | Salinity | Number | ID1 | | | <input type="checkbox"/> |
| 3 | Sample size (g) | Number | Sample size | | | <input type="checkbox"/> |

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The screenshot displays the 'tiamo 2.5 - Alkalinity for 50 ml samples' software interface. On the left, a vertical toolbar contains icons for Workplace, Database, Method, Configuration, and Manual. The main area features a workflow diagram with the following steps: START (Main track), MEAS U (Initial pH-Temperature), CALC (Calcul Voltage pH 4), CALC (Calcul Voltage pH 3), MET U (Titration), and EXPORT (Export to csv). Three overlapping dialog boxes are shown, all titled 'MEAS U - Initial pH-Temperature'. The top-left dialog is in the 'General/Hardware' tab, showing fields for Device name (888_1), Device type (888 Titrande), Sensor (pH electrode), and Stirrer (1). The top-right dialog is in the 'Evaluations' tab, showing checkboxes for 'Fixed endpoint evaluation', 'Minimum evaluation', 'Maximum evaluation', and 'Break point evaluation'. The bottom dialog is in the 'Measuring parameters' tab, showing 'Measurement with drift control' selected, with parameters for Signal drift (1 mV/min), Min. waiting time (0 s), Max. waiting time (30 s), Measuring interval (2.0 s), Stop measured value (off mV), and Temperature (22 °C).

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MET U
Titration

EXPORT
Export to csv

Workplace
Database
Method
Configuration
Manual

CALC - Calcul Voltage pH 4

Command name: Calcul Voltage pH 4

| Result name | Formula | Unit | Decimal places | Assignment | Statistics | Result monitoring |
|----------------|---|------|----------------|------------|-------------------------------------|--------------------------|
| 1 Voltage pH 4 | =(('CV.Voltage TRIS' /1000)-(4- 'CV.pH TRIS' .. | | 2 | RS10 | <input checked="" type="checkbox"/> | <input type="checkbox"/> |

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The screenshot shows the 'tiamo 2.5 - Alkalinity for 50 ml samples' application. The main workflow is as follows:

- START** (Main track)
- MEAS U** (Initial pH-Temperature)
- CALC** (Calcul Voltage pH 4)
- CALC** (Calcul Voltage pH 3)
- MET U** (Titration)
- EXPORT** (Export to csv)


Three dialog boxes are open for configuration:

- CALC - Calcul Voltage pH 3**: Shows a table with columns: Result name, Formula, Unit, Decimal places, Assignment, Statistics, Result monitoring. Row 1: Voltage pH 3, $=('CV.Voltage\ TRIS' / 1000) - (3 - 'CV.pH\ TRIS') * 8,31447215 * ('Initial\ pH-Temperature.ETE' + 273,15) * (ln(10) / 96485,339924) * 1000$, mV, 2, RS09, checked, unchecked.
- Result Voltage pH 3 - RS09**: Shows the formula and unit settings. Unit: mV, Decimal places: 2, Assignment: RS09. Statistics: checked. Description: RS.'Result name'[,VAL] Result value.
- Result Voltage pH 3 - RS09**: Shows options for saving results. Save result as common variable. Common variable: Voltage pH 3. Save result as titer. Save result as global variable. Display message. Record message. Message by e-mail.

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The screenshot shows the 'tiamo 2.5 - Alkalinity for 50 ml samples' application window. The main interface includes a menu bar (File, Edit, View, Insert, Help), a toolbar, and a left sidebar with icons for Workplace, Database, Method, Configuration, and Manual. A central workflow diagram shows the sequence of operations: CALC (Calcul Voltage pH 4), CALC (Calcul Voltage pH 3), MET U (Titration), EXPORT (Export to csv), DATABASE (Transfer to database), and END (Main track). Three configuration dialog boxes for 'MET U - Titration' are open, displaying the following settings:

- Command name:** Titration
- General/Hardware:**
 - Device name: 888_1
 - Device type: 888 Titrando
 - Dosing device: 1
 - Solution: HCl
 - Measuring input: 1
 - Sensor: pH electrode
 - Temperature measurement: automatic
 - Stirrer: 1
 - Stirring rate: 2
 - Switch off automatically
- Start conditions:**
 - Initial measured value:
 - Signal drift: 1 mV/min
 - Min. waiting time: 10 s
 - Max. waiting time: 240 s
 - Start volume:
 - Start volume: 0 mL
 - Dosing rate: maximum mL/min
 - Start measured value:
 - Start measured value: Voltage pH 4' mV
 - Dosing rate: 1 mL/min
 - Start slope:
 - Start slope: off mV/mL
 - Dosing rate: 0.5 mL/min
 - Pause: 120 s
- Titration parameters:**
 - Titration rate: user
 - Measured value acceptance:
 - Signal drift: 1 mV/min
 - Min. waiting time: 0 s
 - Max. waiting time: 154 s
 - Dosing of increments:
 - Volume increment: 0.02 mL
 - Dosing rate: maximum mL/min
 - Temperature: 22 °C

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Command name: Titration

General/Hardware | Start conditions | Titration parameters

Additional evaluations | Additional measured values

Stop conditions | Potentiometric evaluation

Stop conditions

Stop volume: off mL

Stop measured value: Voltage pH 3' mV

Stop EP: off

Volume after EP: off mL

Stop time: off s

Filling rate: maximum mL/min

OK Cancel

MET U - Titration

Command name: Titration

General/Hardware | Start conditions | Titration parameters

Additional evaluations | Additional measured values

Stop conditions | Potentiometric evaluation

Stop conditions

Evaluation without window

EP criterion: 30 mV

EP recognition: all

Evaluation with measured value window (U)

Evaluation with volume window (mL)

OK Cancel

MET U - Titration

Command name: Titration

General/Hardware | Start conditions | Titration parameters | Stop conditions

Potentiometric evaluation | Additional evaluations | **Additional measured values**


Additional calculated measured values

| | Name | Formula | Unit | Decimal places |
|---|-------------|------------------------|------|----------------|
| 1 | Sample size | = 'MV.Sample size (g)' | g | 2 |
| 2 | Salinity | = 'MV.Salinity' | g/L | 2 |

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The screenshot shows the 'tiamo 2.5 - Alkalinity for 50 ml samples' software interface. The main window displays a vertical workflow diagram with the following steps:

- Workplace** (icon: circular arrows)
- Database** (icon: folder)
- Method** (icon: three horizontal bars)
- Configuration** (icon: sliders)
- Manual** (icon: hand holding a coin)

The workflow steps are:


- CALC** (blue box): Calcul Voltage pH 4
- CALC** (blue box): Calcul Voltage pH 3
- MET U** (green box): Titration
- EXPORT** (blue box): Export to csv
- DATABASE** (blue box): Transfer to database
- END** (orange box): Main track

An **EXPORT - Export to csv** dialog box is open, showing the following configuration:

- Command name:
- Export template:
- Buttons:

At the bottom of the software window, it says "Press F1 for help".

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The screenshot shows the 'tiamo 2.5 - Alkalinity for 50 ml samples' software interface. The main window displays a vertical workflow diagram with the following steps:

- Workplace** (Refresh icon)
- Database** (Folder icon)
- Method** (Method icon)
- Configuration** (Sliders icon)
- Manual** (Hand icon)

The workflow steps are:

- CALC** (Calcul Voltage pH 4)
- CALC** (Calcul Voltage pH 3)
- MET U** (Titration)
- EXPORT** (Export to csv)
- DATABASE** (Transfer to database)
- END** (Main track)

A dialog box titled "DATABASE - Transfer to database" is open, showing the following configuration:

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